Gebrauchsanweisung | Instructions for use | Mode d'emploi
Modo de empleo | Modalità d'uso

CE 0483

Fig. 1 Bestandteile der SUS² – Sabbagh Universal Spring.
 Components of the SUS² – Sabbagh Universal Spring.
 Composants du SUS² – Sabbagh Universal Spring.
 Componentes del SUS² – Sabbagh Universal Spring.
 Componenti del SUS² – Sabbagh Universal Spring.

Position für Turbofeder oder Distanzringe
 Position for turbo spring or spacer rings
 Position pour ressort turbo ou anneaux d'écartement
 Posición de el turbo-resorte o anillos distanciadores
 Arresto per molle turbo o anelli distanziatori

Fixierschraube
 Fixing screw
 Vis de fixation
 Tornillo de sujeción
 Vite di fissaggio

Teleskopstange mit Ring
 Telescope rod with ring
 Tige télescopique avec bague
 Varilla telescópica con anillo
 Perno telescopico con anello

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Notizen

Dear Customer

Thank you for choosing this high-quality product from Dentaureum.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at www.dentaureum.de.

1. Manufacturer

Dentaureum GmbH & Co. KG · Turnstraße 31 · 75228 Ispringen · Germany

2. General product description

The SUS² – Sabbagh Universal Spring – produces dentoalveolar adjustments when used as a spring with maximum spring activation (Jasper-Jumper™ effect/elastics substitute).

To correct distoclusion or attain a reposition effect, set the spring activation to minimum (Herbst appliance/activator effect).

Indications: distoclusion/Class II (unilateral or bilateral), molar distalization, diastema closure (aplasia), mandibular joint dysfunction, sleep apnea therapy.

Contraindications: severe protrusion and/or overcrowding in the lower anterior region, poor oral hygiene, gummy smile.

Fig. 2 SUS² in situ.


Fig. 3 Setting the required length or force of the SUS², above as a spring, below as a hinge (Herbst effect).

3. Application of the SUS²

The SUS² as a spring: The SUS² is supplied with the inner spring set at maximum activation (spring path 6.0 mm, force with full compression approx. 2.4 N). This takes into consideration the most common indication of the appliance, i.e. used in the same way as the Jasper-Jumper™ with existing bands/brackets and without laboratory fabricated retention. The effect is comparable to a total of 24 hours of headgear therapy with Class II elastics – mainly a dentoalveolar effect (Fig. 3).

The SUS² as a hinge: The inner spring can be deactivated using the hexagon socket key (REF 607-144-00) by turning the hexagon socket screw of the telescope element as far as it will go, in a counter-clockwise direction (Fig. 3). This converts the flexible SUS² spring as into a rigid SUS² hinge. (Spring path 2.1 mm, buffer force with full compression approx. 1.0 N).

1. Insert the SUS² distally through the eye of the SUS² telescope element (Fig. 4).
2. Insert the SUS² together with the SUS² telescope element distally into the headgear tube and bend it mesially. It is essential to maintain a gap of approx. 4 mm between the head of the and the distal end of the headgear tube (Fig. 5) to ensure that there is enough clearance for the SUS².
3. Bend the SUS² between the band and hook of the buccal tube or shorten it and secure in position using composite/acrylic (**Important!**) (Fig. 6).
4. Bend the SUS² to the distal of the headgear tube by approx. 20° and finely adjust (Fig. 7).
5. Insert the SUS² telescope rod into the SUS² telescope element (guide tube and middle telescope tube) which is fixed in position in the upper jaw.


Fig. 4 Inserting the SUS² into the eye of the SUS² telescope element.

Fig. 5 SUS² in position.

Fig. 6 SUS² bent between the band and hook of buccal tube.

Fig. 7 SUS² bent buccally by 20 to 30 degrees.

Fig. 8 Inserting the SUS² arch adapter on the stainless steel arch.

Fig. 9 L bend.

6. Attach the preassembled SUS² arch adapter (closed side facing upwards) onto the stainless steel archwire from above (minimum 0.41 x 0.56 mm / 16 x 22) between the lower canine and first lower premolar (Fig. 8). The SUS² arch adapter can be also fixed in position between the lower lateral and lower canine if required, e.g. in the case of a missing premolar. L bend (Fig. 9).

Fully tighten the hexagon socket screw using the hexagon socket key (REF 607-144-00). The SUS² arch adapter should not come into contact with the bracket when tightening the screw.

4. Setting and activating the SUS²

One or a combination of the following steps can be taken to ensure that the SUS² meets individual requirements:

The SUS² as a spring: The SUS² turbo spring, with an additional force of approx. 3.0 N, is available if greater activation is required. The turbo spring is placed on the SUS² telescope element so that the three closed coils are pointing towards the SUS² telescope rod (Fig. 10).

SUS² as a hinge: Standard activation is completed using the 1 mm and 2 mm spacer rings supplied (REF 607-103-11 and 607-104-11). Depending on requirements, one or several spacer rings can be clamped on the SUS² telescope rod using special crimping pliers (REF 003-710-00) or Weingart universal pliers (REF 003-120-00) (Fig. 11).

Caution: Only use the SUS² as a hinge on cast or soldered appliances that are retained using bands – the appliance can only be activated using spacer rings (no SUS² turbo spring).


Fig. 10 Additional activation when using as a spring.

Fig. 11 Additional activation when using as a hinge.

5. Tips

Problem	Possible cause	Remedy
Patient bites on the metal or telescope component catches on the premolar brackets	SUS ² not long enough distally, not bent or not adjusted	- SUS² pin should be approx. 4 mm long distally from the headgear tube and bent (see Fig. 5-7) - adjust the SUS² pin buccally or, if required, caudally
Patient cannot bite together properly/SUS ² is too long	SUS ² too long due to extraction or very narrow premolars	- place the SUS² telescope rod/SUS² arch adapter between the lower 2 and lower 3 instead of between the lower 3 and lower 4 - L bend (Fig. 9)
Middle telescope tube jams in the guide tube (no spring function)	(Calculus in the telescope thread) patient with a heavy build-up of calculus	Compensate by using the SUS ² turbo spring, if necessary replace the SUS ² telescope component (Vaseline in the guide tube helps prevent calculus build-up)
The SUS ² telescope rod slides out of the SUS ² telescope element when the mouth is fully open	Patient's mouth opens wider than average, condylar hypermobility	- place the SUS² arch adapter further distally, if possible - explain reinsertion to the patient - physiotherapy/Rotation exercises
Lower anterior teeth are severely protruded	Too little anchorage and/or overactivation	- insert a strong square stainless steel arch min. 0.41 x 0.56 mm / 16 x 22 with lingual crown torque and bend it distally - Elasto-Force plastic chain from 6-6 - temporary anchorage with tomas®-pins (Dentaurum) - reduce the force setting
SUS ² telescope rod catches/rubs during mouth closure	SUS ² telescope rod bent during automatic mouth closure (after yawning)	Replace the SUS ² telescope rod (can be ordered separately)
Brackets loosen frequently	- excessive force (incorrect anchorage/application) - patient bites on the SUS² telescope element (see above)	- rigid hinge and spacer rings (Herbst version/skeletal) should only be anchored on cemented bands similar to Herbst retention - only use the SUS² with activated middle telescope tube on multibracket anchorage (dental)
Fixing screw on the SUS ² arch adapter comes loose	Fixing screw too loose or not tightened enough	Tighten the fixing screw and fix in position using a screw lock adhesive, e.g. Ceka® Bond

6. Composition

The composition is included in the materials list. Refer to www.dentaurum.de.

7. Safety information

The product should not be used, if there is a known allergic reaction to one or more of the material components.

8. Information on single-use products

The SUS² is intended for single use only. Reconditioning a used SUS² (recycling) or its reuse on a patient is not permitted.

9. Availability

SUS² Start-Set Contents:		REF 607-140-01
2 x	SUS ² telescope element	REF 607-147-00
2 x	SUS ² telescope rod, assembled	REF 607-148-00
2 x	SUS ² ball retainer clasps, ø 1.0 mm	REF 607-136-00
2 x	Spacer rings, crimpable, Length 1 mm	REF 607-103-11
2 x	Spacer rings, crimpable, Length 2 mm	REF 607-104-11
2 x	SUS ² turbo spring (Stainless steel), Length 22 mm	REF 607-146-00
1 x	Hexagon socket key	REF 607-144-00

SUS² Set Contents:		REF 607-145-01
2 x	SUS ² telescope element	REF 607-147-00
2 x	SUS ² telescope rod, assembled	REF 607-148-00
2 x	SUS ² ball retainer clasps, ø 1.0 mm	REF 607-136-00
2 x	Spacer rings, crimpable, Length 1 mm	REF 607-103-11
2 x	Spacer rings, crimpable, Length 2 mm	REF 607-104-11
2 x	SUS ² turbo spring (Stainless steel), Length 22 mm	REF 607-146-00

Single Packs		
2 x	SUS ² telescope rod, assembled (each comprising: telescope rod, arch adapter/fixing screw)	REF 607-148-00
2 x	SUS ² telescope rod	REF 607-149-00
2 x	SUS ² arch adapter/fixing screw	REF 607-143-00
1 x	SUS ² telescope element (with hexagon socket screw)	REF 607-147-00
10 x	SUS ² ball retainer clasps, ø 1.0 mm	REF 607-136-00
6 x	SUS ² turbo spring (Stainless steel), Length 22 mm	REF 607-146-00
10 x	Spacer rings, crimpable, Length 1 mm	REF 607-103-11
10 x	Spacer rings, crimpable, Length 2 mm	REF 607-104-11
1 x	Hexagon socket key for Herbst TS / I / II, SUS ²	REF 607-144-00

10. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

11. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at www.dentaurum.de (Explanation of label symbols REF 989-313-00).

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